

PATIENT NAME	Amit Kumar	DATE	18-Jan-2011
REF. BY DR.	Vivek mbbs	SEX / AGE	Male 25 Years
COLLECTION CENTRE	Hospital	LAB NO.	121

RHEUMATOID FACTOR [RF] IN HUMAN SERUM

Result	Positive
Titre	1 : 120
Method	Latex Agalitution
Kit Used	Omega Diagnostics Limited
Controls Used	Positive & Negative

Comments : Examine the test slide under a strong light source after 2 minutes. A positive result is indicated by the obvious agglutination pattern of the latex, in a clear solution. A negative result is indicated by no change in the latex suspension on the test slide.

AVITEX-RF has a detection limit of 8 IU/ml of RF in the patient's serum and the reagent is calibrated against the World Health Organisation (WHO) International Reference Preparation. Positive results will be obtained at a RF serum concentration of 8 IU/ml or more and negative results will be obtained at a RF concentration below 8 IU/ml.